EN SIKKER HÅND

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Referencer:

- 1 Cabana MD, Rand CS, et al Why don't physicians follow clinical practice guidelines? A framework for improvement. *JAMA*. 1999 Oct 20;282(15):1458-65.
- 2 Schuster MA, McGlynn EA, Brook RH. How good is the quality of health care in the United States? *Milbank Q*. 1998;76(4):517-63, 509.
- 3 US Department of Health and Human Services - Agency for Healthcare Research and Quality. Translating Research Into Practice (TRIP)-II Fact Sheet. <http://www.ahrq.gov/research/trip2fac.htm>
- 4 Shimokura G, Weber DJ, et al Factors associated with personal protection equipment use and hand hygiene among hemodialysis staff. *Am J Infect Control*. 2006 Apr;34(3):100-7.
- 5 D De Wandel, D Vogelaers and S Blot. The WHO glove use pyramid: knowledge gaps among Belgian nurses. Presentation at International Conference on Prevention & Infection Control (ICPIC). *BMC Proceedings* 2011, 5(Suppl 6):P114.
- 6 Cutter J, Jordan S. Inter-professional differences in compliance with standard precautions in operating theatres: A multi-site, mixed methods study. *Int J Nurs Stud*. 2012 Mar 27.
- 7 Fuller C, Savage J, Besser S et al "The dirty hand in the latex glove": a study of hand hygiene compliance when gloves are worn. *Infect Control Hosp Epidemiol*. 2011 Dec;32(12):1194-9. Epub 2011 Oct 17.
- 8 Krediet AC, Kalkman CJ, et al Hand-hygiene practices in the operating theatre: an observational study. *Br J Anaesth*. 2011 Oct;107(4):553-8. Epub 2011 Jun 10.
- 9 Girou E, Chai SH, et al Misuse of gloves: the foundation for poor compliance with hand hygiene and potential for microbial transmission? *J Hosp Infect*. 2004 Jun;57(2):162-9.
- 10 Parmeggiani C, Abbate R et al Healthcare workers and health care-associated infections: knowledge, attitudes, and behavior in emergency departments in Italy. *BMC Infect Dis*. 2010 Feb 23;10:35.
- 11 Stein AD, Makarawo TP, Ahmad MF. A survey of doctors' and nurses' knowledge, attitudes and compliance with infection control guidelines in Birmingham teaching hospitals. *J Hosp Infect*. 2003 May;54(1):68-73.
- 12 Eveillard M, Raymond F, et al Impact of a multi-faceted training intervention on the improvement of hand hygiene and gloving practices in four healthcare settings including nursing homes, acute-care geriatric wards and physical rehabilitation units. *J Clin Nurs*. 2011 Oct;20(19-20):2744-51. doi: 10.1111/j.1365-2702.2011.03704.x. Epub 2011 Mar 3.
- 13 Howard DP, Williams C, et al A simple effective clean practice protocol significantly improves hand decontamination and infection control measures in the acute surgical setting. *Infection*. 2009 Feb;37(1):34-8. Epub 2008 Dec 5.
- 14 Trick WE, Vernon MD, et al Multicenter intervention program to increase adherence to hand hygiene recommendations and glove use and to reduce the incidence of antimicrobial resistance. *Infect Control Hosp Epidemiol*. 2007 Jan;28(1):42-9. Epub 2007 Jan 3.
- 15 Twomey CL. Double gloving: a risk reduction strategy. *Jt Comm J Qual Saf*. 2003 Jul;29(7):369-78.
- 16 <http://www.aaos.org/about/papers/advistmt/1018.asp>
- 17 http://www.cdc.gov/sharpsafety/pdf/sharpsworkbook_2008.pdf
- 18 Statement on Sharps safety [2007] http://www.facs.org/fellows_info/statements/st-58.html
- 19 http://www.aorn.org/Books_and_Publications/Perioperative_Standards_and_Recommended_Practices/Perioperative_Standards_and_Recommended_Practices.aspx#axzz20S20hrKH
- 20 Tanner J, Parkinson H. Double gloving to reduce surgical cross-infection. *Cochrane Database Syst Rev*. 2006 Jul 19;(3):CD003087
- 21 Guo YP, Wong PM, Li Y, Or PP. Is double-gloving really protective? A comparison between the glove perforation rate among perioperative nurses with single and double gloves during surgery. *Am J Surg*. 2012 Feb 16.
- 22 Korniewicz D, El-Masri M. Exploring the benefits of double gloving during surgery. *AORN J*. 2012 Mar;95(3):328-36. Association of Perioperative Registered Nurses (AORN)
- 23 Salimone LA. Mandating double gloving. *AORN Journal Volume 95, Issue 3, Pages 328-336, March 2012.*
- 24 Twomey CL. Double gloving: a risk reduction strategy. *Jt Comm J Qual Saf*. 2003 Jul;29(7):369-78.
- 25 Sackett DL, Rosenberg W et al Evidence based medicine: what it is and what it isn't. *British Medical Journal* 1996;312:71

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
Referencer, se sidste side.

Referencer:

1. Tanner J, Parkinson H. Double gloving to reduce surgical cross-infection. *Cochrane Database of Systematic Reviews* 2006, Issue 3. Art. No.: CD003087. DOI: 10.1002/14651858.CD003087.pub2. Review content assessed as up-to-date: 28 July 2009.
2. Hollaus P. Et al. Glove perforation rate in open lung surgery. *European Journal of Cardio-thoracic Surgery* 15;(1999):461-464.
3. Hübner NO. Et al. Migration through punctured surgical gloves under real surgical conditions *BMC Infectious Diseases* 2010;10:192-197.
4. Esterban JI. Et al. Transmission of hepatitis C virus by a cardiac surgeon. *The new England Journal of Medicine* 1996;334(9):555-560.
5. Harpaz R. Et al. Transmission of hepatitis B virus to multiple patients from a surgeon without evidence of inadequate infection control. *N Engl J Med* 1996;334:549-54.
6. McNeil SA. Et al. Outbreak of sternal surgical site infections due to *Pseudomonas aeruginosa* traced to a scrub nurse with onychomycosis. *Clin Infect Dis Epub* 2001;33(3):317-323.
7. Van Den Broek PJ. Epidemic of prosthetic valve endocarditis caused by *Staphylococcus epidermidis*. *British Medical Journal* 1985;291:949-950.
8. Wooster DL. Et al. Intraoperative bacterial contamination of vascular grafts: a prospective study. *Can J Surg* 1985; 28(5):407-409.
9. Olsen K. Et al. Increased risk of transmission of hepatitis C in open heart surgery compared with vascular and pulmonary surgery. *Ann Thorac Surg*. 2010;90(5):1425-1431.
10. Bennett NT, Howard RJ. Quantity of blood inoculated in a needlestick injury from suture needles. *Journal of the American College of Surgeons* 1994; 178(2):107-110.
11. Preventing the Transmission of Bloodborne Pathogens. Revised June 2008 American Academy of Orthopaedic Surgeons (AAOS) page 5. (Double gloving is recommended)
12. Workbook for Designing, Implementing and Evaluating a Sharps Injury Prevention Program. Centers for Disease Control and Prevention 2008. (p.19 Work-practice control in the operating room = double gloving)
13. Statement on Sharps safety. American College of Surgeons 2007. (The ACS recommends the universal adoption of the double glove technique in order to reduce body fluid exposure caused by glove tears and sharps injuries in surgeons and scrub personnel)
14. Perioperative Standards and Recommended Practices. AORN 2011 Recommenden practice for prevention of transmissible infections in the perioperative practice setting. (p. 293, Recommendation IV. Health care practitioners should double-glove during invasive procedures).
15. Perioperative Standards and Recommended Practices. AORN 2011 Sharp injuries prevention in the perioperative setting. (p. 640, Perioperative-specific risk reduction strategies = During the operative or other invasive procedure: wear two pairs of gloves [ie. Double gloving]).
16. Korniewicz D, El-Masri M. Exploring the benefits of double gloving during surgery. *AORN J*. 2012 Mar;95(3):328-36.

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